

Abbott – Recall of FreeStyle Libre® 3 and FreeStyle Libre 3 Plus Sensors

- On November 24, 2025, Abbott announced a voluntary recall of multiple lots of FreeStyle Libre 3 and FreeStyle Libre 3 Plus sensors because the sensors may provide incorrect low glucose readings.
- The product was distributed nationwide.

Product Description	NDC number	Lot number
FreeStyle Libre 3 sensor	57599-0818-00	T60003054; T60003092; T60003159 T60003374; T60003434; T60003547 T60003085; T60003099; T60003160 T60003375; T60003507; T60003564 T60003088; T60003113; T60003271 T60003391; T60003534; T60003628 T60003089; T60003136; T60003284 T60003426; T60003546; T60003646
FreeStyle Libre 3 Plus sensor	57599-0844-00 57599-0844-01	T60002984; T60003366; T60003455; T60003492; T60003560; T60003592; T60002994; T60003369; T60003460; T60003494; T60003562; T60003596; T60002995; T60003377; T60003470; T60003496; T60003563; T60003601; T60002996; T60003378; T60003471; T60003497; T60003567; T60003611; T60003007; T60003379; T60003473; T60003498; T60003568; T60003617; T60003009; T60003384; T60003478; T60003518; T60003570; T60003618; T60003010; T60003398; T60003479; T60003542; T60003573; T60003619; T60003138; T60003400; T60003480; T60003543; T60003575; T60003631; T60003243; T60003404; T60003481; T60003544; T60003576; T60003639; T60003255; T60003446; T60003488; T60003550; T60003585; T60003641; T60003348; T60003449; T60003490; T60003552; T60003589; T60003683; T60003352; T60003454; T60003491

- FreeStyle Libre 3 and FreeStyle Libre 3 Plus sensors provide continuous glucose monitoring levels for patients with diabetes. The FreeStyle Libre 3 system uses a special sensor that reads blood sugar (glucose) levels every minute and sends the measurements to a smart phone app. The readings tell people with diabetes how to

manage their condition – it may be in the correct range (no action needed), too high (insulin may be needed), or too low (eating food may be needed).

- The FreeStyle Libre 3 and FreeStyle Libre 3 Plus sensor is attached to the back of the upper arm with a special disposable applicator, a small filament is inserted under the skin, and an adhesive pad holds it in place.
- Using the recalled sensors may lead to incorrect low glucose readings. Incorrect low glucose readings over an extended period may lead to wrong treatment decisions. These decisions may pose serious health risks, including potential injury or death, or other less serious complications.

Please contact your doctor immediately if you experience any problems that you think might be related to using the recalled sensors.

- Visit www.FreeStyleCheck.com to confirm if your sensor is impacted and to receive a free replacement from Abbott.

For more information about this recall, please contact your pharmacy or the following company:

Company	Contact information
Abbott	1-833-815-4273 (8:00 a.m. – 8:00 p.m. ET, seven days a week).